Analytical Validation of the Signatera[™] RUO Assay, a Highly Sensitive Patient-Specific Multiplex PCR NGS-Based Non-Invasive Cancer Recurrence Detection and Therapy Monitoring Assay

Himanshu Sethi, Raheleh Salari, Ryan Swenerton, Svetlana Shchegrova, Dina Hafez, Prashanthi Natarajan, Scott Dashner, Tony Tin, Mustafa Balcioglu, Hsin-Ta Wu, Styrmir Sigurjonsson, Bernhard Zimmermann Natera Inc., San Carlos, CA, USA

Introduction

- The identification of tumor mutations in circulating cell-free DNA (cfDNA) holds great potential for the noninvasive detection of cancer relapse before clinical manifestation, detection of minimal residual disease after curative-intent treatment, and detection of therapeutically relevant mutations.
- Previously, we demonstrated 93% sensitivity, with no false positives, for early relapse detection in patients with early stage lung cancer using a proof-of-concept version of the Signatera[™] assay.¹

Objective

• To review and report analytical validation results for the detection of tumor-specific variants by the current version of the assay

Methods


Signatera RUO

• The Signatera (RUO) process starts with identifying and prioritizing somatic mutations from whole-exome sequencing of tumor and matched normal samples. Patient-specific multiplex-PCR assays targeting 16 somatic single-nucleotide and indel variants are then assayed by massively parallel sequencing in plasma samples collected throughout the patient's disease course to help detect and monitor circulating tumor DNA (**Figure 1**).

Analytical Validation

- The analytical validation for the current version of the Signatera (RUO) assay was performed on two breast cancer cell lines (HCC2218, HCC1395), one lung cancer cell line (NCI-H1395) and their matched normal counterparts (HCC2218-BL, HCC1395-BL, and NCI-H1395-BL, respectively). Varying amounts of tumor cell line DNA (0%, 0.005%, 0.01%, 0.03%, 0.05%, 0.1%, 0.3%, 0.5%, 1%) were titrated into its respective matched normal cell line DNA.
- Multiplex-PCR assay primer pools (each consisting of 16 primer-pair assays specific to high-confidence somatic mutations) were designed using whole exome data from the corresponding tumor cell line DNA and its matched normal cell line DNA.
- The starting total input into library prep for each reaction was 20K genome equivalents; SNV and indel targets from the corresponding tumor DNA spike-in samples were amplified using the multiplex-PCR assay primer pools mentioned above.
- The mPCR products were barcoded, then pooled with other mPCR barcoded products, and subsequently sequenced on an Illumina HiSeq 2500 Rapid Run with 50 cycles of paired-end reads using the Illumina Paired End v2 kit with an average read depth of ~100,000/assay.

Figure 1. Signatera (RUO) Process




DAYS SINCE SURGERY

Steps 1 and 2 are optional if the whole-exome sequencing data are available.

Results

• Figure 2 depicts sequencing QC process checks across 4 HiSeq PE 2x50 runs, including, background transition and transversion error rates, and average per target depth-of-reads (DOR) of ~100K (where any targets receiving less than 5K reads were considered failed and not taken into account for calling).





• An analytical sensitivity of ~67% for SNV detection at 0.03% spiked-in tumor DNA was achieved with Signatera (RUO) (**Table 1**).

Table 1. Analytical Sensitivity Results for Signatera (RUO)					
Tumor DNA Concentration (%)	Average Mutant Copies*	Detected SNVs	Total SNVs Used	Sensitivity (%)	95% CI
0.005	1	120	794	15.1	12.7 - 17.8
0.010	2	246	635	38.7	34.9 - 42.7
0.030	6	291	433	67.2	62.6 - 71.6
0.050	10	208	244	85.3	80.2 - 89.5
0.100	20	227	228	99.6	97.6 - 100
0.300	60	177	177	100	98 - 100
0.500	100	120	120	100	97 - 100
1.000	200	76	76	100	93 - 100
*Based on an average number of 20K starting genome copies, as expected from 10 ml of plasma.					

- Estimated sample-level sensitivities for Signatera (RUO) when at least two SNVs are called from a set of • For a given set of SNVs, the number of expected input percentages versus detected variant allele frequency (VAF) percentages for the differing spiked-in mutant DNA concentrations is shown across six targets, displaying 16 high-confidence SNVs, are listed in Table 2. high sensitivity above 0.03% tumor DNA concentrations (Figure 4).
- The false negatives shown at ≤0.01% SNV input (starting with <2 mutant copies) include mutant molecules lost due to sampling (**Figure 4**).

Abstract Number 4542

Table 2. Estimated Sample-Level Sensitivity				
Tumor DNA Concentration (%)	Sensitivity (%)			
0.005	49.97 -72.02			
0.010	94.55 - 99.56			
0.030	99.97 - 99.99			
0.050	100			
0.100	100			
0.300	100			
0.500	100			
1.000	100			

Conclusions

- The Signatera RUO assay provides a novel method of non-invasively detecting recurrence of personalized cancer signatures in plasma by ultra-deep sequencing of custom-made multiplex PCR assays (selected from a patient's tumor) up to 14.7 months in colorectal cancer (see Abstract Number 1590) with high sensitivity, high specificity, and low error rates.
- Based on the analytical validation results, the Signatera RUO assay, on an SNVlevel, has a 99.9% specificity and a greater than 65% sensitivity above 0.03% tumor fraction and a 100% sensitivity above 0.1% tumor fraction.
- On a sample-level, the Signatera RUO assay has greater than 95% sensitivity at 0.01% tumor fraction, nearly 100% sensitivity at 0.03% tumor fraction, and 100% sensitivity at 0.05% and above tumor fraction.
- These data demonstrate high rate of detection at single molecule mutant levels; they also suggest that lower plasma volumes may be utilized to achieve the same single-molecule detection with high specificity.
- The performance of the Signatera assay suggests the potential for it to determine chemotherapy treatment effectiveness.





