

Highlights from Prospera's clinical validation

Prospera™ is clinically validated to assess active rejection with great accuracy in the largest published donor-derived cell-free DNA (dd-cfDNA) validation study on renal transplant patients.

Published in *Journal of Clinical Medicine*, 2018



COHORT

Diverse patient population studied with a trusted partner

- Partnered with the University of California, San Francisco
- 217 biopsy-matched renal allograft samples
- Variety of ethnic and racial demographics
 - Hispanic/Latino (n=50)
 - Caucasian (n=74)
 - African American (n=31)
 - Asian (n=31)



DESIGN

Blinded, retrospective analysis using gold standard of kidney biopsy as truth

- Blinded pathology readings and analysis of blood samples
- Included same-day biopsy matched samples in primary analysis
- Reviewed by a single pathologist to ensure consistency and quality in the data
- Used prospectively selected, previously published cut-off of 1% dd-cfDNA to discriminate between active rejection and not active rejection




ANALYSIS

Classification based on recent guidelines

- Assessed active rejection in the following categories: antibody-mediated rejection (ABMR), T cell-mediated rejection (TCMR \geq 1a) and mixed (ABMR/TCMR)
- Categorized biopsies with “borderline rejection” in the non-rejection group instead of active rejection, consistent with the latest BANFF 2017 criteria

Clinical Validation Results

More sensitive and specific than standard screening tools¹




Subanalysis:

Consistent performance regardless age

- Similar performance of sensitivity (89%), specificity (76%) and AUC (0.88) an **adults-only cohort** to the overall study population⁴
- Significantly higher predictive power (AUC of 0.84) compared to serum creatinine (AUC of 0.58) in assessing active rejection in an **adults-only, for-cause sub-cohort**. This serum creatinine performance of 0.58 AUC is in-line with other previously published serum creatinine performance metrics (0.542², 0.633³)

Robust to all rejection types



- In contrast to other dd-cfDNA studies, Natera was able to differentiate T cell-mediated rejection from no rejection
- Median dd-cfDNA did not differ significantly between ABMR, TCMR and mixed rejection

Similar performance in stable and those with clinical signs of rejection

	Sensitivity	Negative predictive value*
Clinical AR	86%	93.8%
Subclinical AR	92%	98.3%**

* Assuming 25% AR Prevalence (higher risk population)

** Assuming 10% AR Prevalence in a low risk population

- Prospera is the first commercial cfDNA assay to publish performance in renal transplant surveillance situations, providing results that can enable physicians to manage patients with previously unsuspected rejection^{1,2}

References:

1. Sigdel TK, Archila FA, Constantin T, et al. Optimizing detection of kidney transplant injury by assessment of donor-derived cell-free DNA via massively multiplex PCR. *Journal of Clinical Medicine*. 2019;8(1):19.
2. Bloom RD, et al. *J Am Soc Nephrol*. 2017 Jul;28(7):2221-2232.
3. Kaplan, et al. *Am Journal of Transplantation*. 2003; 3:1560-1565
4. Sigdel TK, Acosta Archila F, Navarro S, et al. Rapid Detection of Kidney Transplant Injury by Quantifying Donor-Derived Cell-Free DNA via Massively Multiplex PCR. Poster presented at: American Transplant Congress (ATC); 2019 June 01-05, Boston, MA.

201 Industrial Road, Suite 410 | San Carlos, CA 94070 | www.natera.com | 1.650.249.9090 | Fax 1.650.730.2272

This test was developed by Natera, Inc., a laboratory certified under the Clinical Laboratory Improvement Amendments (CLIA). This test has not been cleared or approved by the US Food and Drug Administration (FDA). Although FDA does not currently clear or approve laboratory-developed tests in the US, certification of the laboratory is required under CLIA to ensure the quality and validity of the tests.

© 2019 Natera, Inc. All Rights Reserved. PRO_OS_SigdelCVHighlights_20190911_NAT-801996

 **natera**[®]
Conceive. Deliver. Thrive.