



Become a Leader in Genetic Testing

Tap into the Power of Natera Technology

Natera's Constellation Platform

Constellation™ allows laboratories worldwide to access key technologies that Natera uses in its own lab. Partner labs now have the opportunity to develop, validate, and clinically offer complex genomics assays in their own facilities.



Sample
Collection



DNA
Extraction



SNP Probe &
Library Preparation



Sequencing



Data Analyst
& Test Report

Present

Panorama™ Non-Invasive
Prenatal Test

Future

Liquid Biopsies for
Oncology Applications

Natera is the Market Leader in NIPT



- Highly accurate screening for detection of aneuploidy, triploid, molar and vanished twin pregnancies, and specific microdeletions
- Accurate fetal sex determination
- Significant fewer false positives than Massively Parallel Shotgun Sequencing (MPSS)
- Fetal fraction analysis

Panorama Clinical Performance

	Sensitivity			Specificity		
Trisomy 21	83/83	>99%	(CI: 95.6–100%)	1,108/1,108	>99%	(CI: 99.7–100%)
Trisomy 18	27/28	96.4%	(CI: 81.7–99.9%)	1,164/1,165	>99%	(CI: 99.5–100%)
Trisomy 13	13/13	>99%	(CI: 75.3–100%)	1,180/1,180	>99%	(CI: 99.7–100%)
Monosomy X	13/14	92.9%	(CI: 66.1–99.8%)	1,179/1,180	>99%	(CI: 99.5–100%)
Triploidy	8/8	99%	(CI: 63.1–100%)	272/272	>99%	(CI: 98.7–100%)
Female	469/469	>99%	(CI: 99.2–100%)	533/533	>99%	(CI: 99.3–100%)
Male	533/533	>99%	(CI: 99.3–100%)	469/469	>99%	(CI: 99.2–100%)

Natera's Constellation Offering to Partner Labs



NIPT Test Package

- 1. CE Marked Reagents:**
Library Prep and PCR Multiplex Kits
- 2. CE Marked Bioinformatics:**
Streamlined and optimized for Panorama
- 3. Standard Operating Protocols:**
Protocols based on best practices employed by Natera's CLIA lab



Not Available in the USA

The Natera Advantage



Cost efficient and clinically relevant targeted assays



Panorama™ protocol easy to implement




Partner can access same algorithm in their own labs



Long-term and expanding strategic partnerships in clinical genomics

- Dedicated scientific staff to support site preparation and workflow customization
- Step-by-step implementation program with live training
- Assistance in procurement of validation samples
- Dedicated software integration team
- Certification program to recognize laboratory, operator, and process quality
- Field-based technical support team for on-market support
- Natera branding and global marketing programs



For further information about Constellation

Web: www.Natera.com

Email: ApplicationSupport@natera.com

Not Available in the USA



201 Industrial Road, Suite 410 | San Carlos, CA 94070 | 1-650-249-9090 | Fax 1-650-730-2272



This test was developed by Natera, Inc., a laboratory certified under the Clinical Laboratory Improvement Amendments (CLIA). This test has not been cleared or approved by the U.S. Food and Drug Administration (FDA). Although FDA does not currently clear or approve laboratory-developed tests in the U.S., certification of the laboratory test is required under CLIA to ensure the quality and validity of the tests.

CAP accredited and CLIA certified, ISO 13485. © 2016 Natera, Inc. All Rights Reserved. CONSTELLATION-MD (5-2016) INTL NAT-801031